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
Fusion Tags Available for pET Constructs

Tag	N/C Terminal or Internal (I)	Size (aa)	Basis for Detection and/or Purification	Applications	pET Vector Series
Dsb•Tag	N	208 (DsbA) 236 (DsbC)	potential periplasmic localization	soluble protein, periplasmic disulfide bond formation, isomerization	39, 40
GST•Tag	N	220	glutathione affinity, mAb, enzymatic activity	purification, Western blot, quantitative assay	41, 42, 49, 60-DEST
His•Tag	N, C, I	6, 8, or 10	mAb, metal chelation chromatography (native or denaturing)	His•Bind® resin purification, Western blot	14-16, 19-52, (53-62)-DEST
HSV•Tag	C	11	mAb	Western blot, immunofluorescence	25, 27, 43.1, 44
KSI	N	125	highly expressed hydrophobic domain	small protein/peptide production/purification, insoluble protein	31
Nus•Tag	N	495	mAb	soluble protein, cytoplasmic disulfide bond formation in <i>trxB</i> ⁻ hosts, Western blot	43.1, 44, 50, (57, 58)-DEST
<i>pelB</i>	N	20	potential periplasmic localization	protein export/folding	20, 22, 25, 26, 27
S•Tag	N, I	15	mAb, S-protein (104 aa) affinity	Western blot, quantitative assay, purification	29, 30, 32, 39, 40, 41, 42, 43.1, 47, 48, 49, 50, (54, 56, 58)-DEST
Strep•Tag II	N	8	mAb, Strep•Tactin affinity	purification, Western blot	51, 52, (53, 55, 57, 59, 61, 62)-DEST
T7•Tag	N, I	11	mAb	Western blot, immunoprecipitation, purification	3, 9, 11, 17, 21, 23, 24, 28
Trx•Tag	N	109	mAb, thioredoxin	soluble protein, cytoplasmic disulfide bond formation in <i>trxB</i> ⁻ and <i>gor</i> ⁻ hosts, Western blot	32, 48, 59-DEST

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





Analysis of pET-32a(+) IL-2 Clones Induced at 15°, 25° and 37°

A. Soluble fraction

B. Insoluble fraction

University of Oklahoma
School of Chemical Engineering and Materials Science
Recombinant Protein Solubility Prediction

Type (or cut and paste) your protein sequence below, click on the "Submit" button, and the solubility probability of your protein will be calculated. The statistical model predicts protein solubility **assuming the protein is being overexpressed in *Escherichia coli***. If there are numbers, spaces, or other characters in your sequence, don't worry, they won't affect the calculation. For more information on the solubility model used here, see the [references](#) below.

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




* Copy number estimates are based on agarose gel analysis (3).

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Veggie™ BL21(DE3)
Veggie BL21(DE3)pLysS

Insoluble protein/ No activity

Reason: Reduction of disulfide bonds resulting in misfolded protein

Origami™ 2
Origami 2(DE3)
Rosetta-gami™ 2
Rosetta-gami 2(DE3)
Rosetta-gami B
Rosetta-gami B(DE3)

Solution:
Minimize protein reduction in cytoplasm; use *trxB/gor* hosts

Reason: High levels of expression resulting in misfolded protein

Tuner™
Tuner(DE3)
Rosetta-gami B
Rosetta-gami B(DE3)

Solution:
Attenuate expression/titrate IPTG; use *LacY*⁻ hosts

Truncated protein

Reason: *E. coli* codon bias

Rosetta™
Rosetta(DE3)
Rosetta 2
Rosetta 2(DE3)
Rosetta-gami™ 2
Rosetta-gami 2(DE3)
Rosetta-gami B
Rosetta-gami B(DE3)
RosettaBlue™
RosettaBlue(DE3)

Solution:
Use a host that supplies rare tRNAs

Stabilizing target plasmids

Reason: Target plasmid unstable do to repetitive sequences

BLR(DE3)
HMS174
HMS174(DE3)
NovaBlue
NovaBlue(DE3)

Solution:
Use *recA*⁻ hosts

Toxic protein

Symptom: No protein/cell death

Tuner™
Tuner(DE3)
NovaBlue
NovaBlue(DE3)
Any pLysS host

Solution:
Suppress basal expression; use a stringent control host

Protein labeling

B834
B834(DE3)

Use a methionine auxotroph host for higher specific activity

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



Tissue plasminogen activator (9 disulfide bonds) was expressed in three different hosts. 10 μ g of soluble protein was spotted onto fibrin agarose plates and incubated for 24 h at 37° C. Zone of clearing measures biological activity of tPA.

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
Codon usage is expressed as the fraction of all possible codons for a given amino acid. "All genes" is the fraction represented in all 4,290 coding sequences in the *E. coli* genome (6). "Class II" is the fraction represented in 195 genes highly and continuously expressed during exponential growth (7).

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
← β-gal (5 × CGG mutant)

Lane	Sample
1	Rosetta(DE3) total cell extract (induced)
2	Rosetta 2(DE3) total cell extract (induced)
M	Perfect Protein™ Markers, 10–225 kDa

Effect of consecutive CGG rare codons on target protein expression

A pET-15b recombinant plasmid containing five consecutive CGG codons near the 5' end of the β-gal coding region was transformed into Rosetta™(DE3) and Rosetta 2(DE3) cells. The cells were induced with IPTG for 3 hr and harvested by centrifugation.

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Strain	DrugR	Genotype	Description
Rosetta™ ²	cam	F- ompT hsdSB(rB-mB-) gal dcm lacY1 pRARE	protease deficient, lac permease deficient, rare tRNAs
RosettaBlue™	cam, tet	endA1 hsdR17 (rK12-mK12 +) supE44 thi-1 recA1 gyrA96 relA1 lac F'[proA+B+ lacIqZDM15::Tn10 pRARE	K12, rare tRNAs
Origami™ ²	tet, Str, kan	Dara-leu7697 DlacX74 DphoAPvu II phoR araD139 ahpC galE galK rpsL F'[lac+(lacIq)pro] gor522::Tn10 trxB::Str,Tet	K12, oxidizing cytoplasm
Origami B	tet, kan	F- ompT hsdSB(rB-mB-) gal dcm lacY1 ahpC gor522::Tn10 trxB::kan	protease deficient, lac permease deficient, oxidizing cytoplasm
Rosetta-gami™	tet, kan, cam	<i>Δara-leu7697 ΔlacX74 ΔphoAPvu II phoR araD139 ahpC galE galK rpsL F'[lac⁺(lacI^q)pro] gor522::Tn 10 trxB::kan pRARE</i>	K12, rare tRNAs, oxidizing cytoplasm
Rosetta-gami B	tet, kan, cam	F ⁻ ompT hsdS _B (r _B ⁻ m _B ⁻) gal dcm lacY1 ahpC gor522::Tn 10 trxB::kan pRARE	protease deficient, lac permease deficient, rare tRNAs, oxidizing cytoplasm

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Veggie BL21(DE3) and Veggie BL21(DE3)pLysS

rLysozyme Solution, Veggie™ Grade

Animal-Free Protease Inhibitor Cocktails

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